

Quantitative Laser Diagnostics for Combustion Chemistry and Propulsion

Lecture 12: Shock Tube Applications with Lasers

1. Laser Absorption Theory
2. Survey of Capabilities
3. Kinetics Applications:
 - Rate Constant Measurements
 - Multi-Species Time-Histories
4. New Species Diagnostics


1. Laser Absorption Theory

- Governing Equation: Beer-Lambert Law


$$I_v/I_0 = \exp(-S_{lu} \Phi(\nu) X_{\text{species}} P_{\text{total}} L)$$

Line
Strength

Line
shape



Lineshape Function: $\Phi(\nu)$



- Quantitative absorption requires database for S , Φ
- What species have been measured?

2. Survey of Capabilities: Species and Wavelengths

Ultraviolet

CH ₃	216 nm
NO	225 nm
O ₂	227 nm
HO ₂	230 nm
OH/CH ₂ O	306 nm
NH	336 nm

Visible


CN	388 nm
CH	431 nm
NCO	440 nm
NO ₂	472 nm
NH ₂	597 nm
HCO	614 nm

Infrared

CO	2.3 μm
H ₂ O	2.5 μm
CO ₂	2.7 μm
Fuel/CH ₂ O	3.4 μm
NO	5.2 μm
C ₂ H ₄	10.5 μm
C ₃ H ₆ /i-C ₄ H ₈	11/11.3 μm

**First use of tunable dye lasers in
shock tubes (1982)**

Spectra-Physics 380



2. Survey of Capabilities: Species and Wavelengths

Ultraviolet

CH ₃	216 nm
NO	225 nm
O ₂	227 nm
HO ₂	230 nm
OH/CH ₂ O	306 nm
NH	336 nm

Visible


CN	388 nm
CH	431 nm
NCO	440 nm
NO ₂	472 nm
NH ₂	597 nm
HCO	614 nm

Infrared

CO	2.3 μm
H ₂ O	2.5 μm
CO ₂	2.7 μm
Fuel/CH ₂ O	3.4 μm
NO	5.2 μm
C ₂ H ₄	10.5 μm
C ₃ H ₆ /i-C ₄ H ₈	11/11.3 μm

Ultra-fast lasers used to extend UV tuning range (2009)

Coherent MIRA
Ti-Sapphire Ring Laser



2. Survey of Capabilities: Species and Wavelengths

Ultraviolet

CH ₃	216 nm
NO	225 nm
O ₂	227 nm
HO ₂	230 nm
OH/CH ₂ O	306 nm
NH	336 nm

Visible

CN	388 nm
CH	431 nm
NCO	440 nm
NO ₂	472 nm
NH ₂	597 nm
HCO	614 nm

Infrared

CO	2.3 μm
H ₂ O	2.5 μm
CO ₂	2.7 μm
Fuel/CH ₂ O	3.4 μm
NO	5.2 μm
C ₂ H ₄	10.5 μm
C ₃ H ₆ /i-C ₄ H ₈	11/11.3 μm

**New lasers allow simple access to
mid-IR**

M² UV Ti:Sapp Laser



2. Survey of Capabilities: Species and Wavelengths

Ultraviolet

CH ₃	216 nm
NO	225 nm
O ₂	227 nm
HO ₂	230 nm
OH/CH ₂ O	306 nm
NH	336 nm

Visible

CN	388 nm
CH	431 nm
NCO	440 nm
NO ₂	472 nm
NH ₂	597 nm
HCO	614 nm


Infrared

CO	2.3 μm
H ₂ O	2.5 μm
CO ₂	2.7 μm
Fuel/CH ₂ O	3.4 μm
NO	5.2 μm
C ₂ H ₄	10.5 μm
C ₃ H ₆ /i-C ₄ H ₈	11/11.3 μm

**How sensitive are laser absorption
diagnostics in shock tubes?**

2. Laser Absorption Yields High Sensitivity


Representative Detection Limits: Large Molecules



- Large molecules: 10-100's ppm

2. Laser Absorption Yields High Sensitivity


Representative Detection Limits: Diatomic Molecules



- **Diatomic molecules @ 1500K:**
 - sub-ppm detectivity for UV absorbers
 - ppm detectivity for IR absorbers

3. Kinetics Applications

1. Foundation Fuel Kinetics
 $\text{OH} + \text{H}_2 = \text{H}_2\text{O} + \text{H}$
2. Methyl Ester Kinetics
 $\text{ME} + \text{OH} = \text{Products}$
Methyl Formate pyrolysis
3. Butanol Isomer Kinetics
n-Butanol pyrolysis
t-butanol + OH using isotopic labeling



3.1 Foundation Fuel Kinetics:

- **Recent Elementary Rate Constant Measurements using Shock Tube/Laser Absorption Methods**

- $\text{H} + \text{O}_2 = \text{OH} + \text{O}$
- $\text{H}_2\text{O}_2 + \text{M} = \text{OH} + \text{OH} + \text{M}$
- $\text{OH} + \text{H}_2\text{O}_2 = \text{H}_2\text{O} + \text{HO}_2$
- $\text{OH} + \text{HO}_2 = \text{H}_2\text{O} + \text{O}_2$
- $\text{HO}_2 + \text{HO}_2 = \text{H}_2\text{O}_2 + \text{O}_2$
- $\text{CH}_3 + \text{HO}_2 = \text{Products}$
- $\text{C}_2\text{H}_4 + \text{M} = \text{C}_2\text{H}_2 + \text{H}_2 + \text{M}$
- $\text{OH} + \text{H}_2 = \text{H}_2\text{O} + \text{H}$




3.1 Foundation Fuel Kinetics

How Do We Measure Individual Reaction Rates?

- Two-Step Process

1. Design: use sensitivity analysis to kinetically isolate reactions
2. Execute: use shock tubes for step heating and laser absorption for species detection




- Accurate k values critical to combustion modeling

3.1 Foundation Fuel Kinetics


- **Example: $\text{H} + \text{O}_2 \rightarrow \text{OH} + \text{O}$**
Rate Measurement using H_2O Laser at 2.55 mm

Shock Conditions: 1472 K, 1.8 atm, 0.1% O_2 /0.9% H_2 /Ar

H_2O Sensitivity




H_2O Time-History



- **H_2O time-history provides precise determination of k_1**


3.1 Foundation Fuel Kinetics

Arrhenius Plot: $\text{H} + \text{O}_2 \rightarrow \text{OH} + \text{O}$



3.1 Foundation Fuel Kinetics

Arrhenius Plot: $\text{H} + \text{O}_2 \rightarrow \text{OH} + \text{O}$



□ Modern laser methods significantly reduce measurement uncertainty!

Next example: $\text{OH} + \text{H}_2 \rightarrow \text{H}_2\text{O} + \text{H}$

3.1 Foundation Fuel Kinetics


- Elementary Reaction Rate Determination:



Motivation: Large uncertainty in $k_{\text{OH}+\text{H}_2}$ gives large uncertainty in modeled H_2 /air flame speeds

Experimental Strategy: Direct determination of rate constant

- Pseudo-first order experiment
- Fuel in excess
- TBHP used as a prompt OH precursor
 - Useful T range (850 to 1350 K)
 - Pioneered by Bott and Cohen (1984), also used at Argonne




TBHP: tert-butylhydroperoxide

Acetone

3.1 Foundation Fuel Kinetics


- Representative OH Absorption Data



- High SNR data, ppm sensitivity
- $k_{\text{best-fit}}$ determined within 3-5%

3.1 Foundation Fuel Kinetics


- Arrhenius Plot: $\text{OH} + \text{H}_2 \rightarrow \text{H}_2\text{O} + \text{H}$



- Very low overall uncertainty in k : $\pm 17\%$ (2σ)

3.1 Foundation Fuel Kinetics

- Comparison with Past Work: $\text{OH} + \text{H}_2 \rightarrow \text{H}_2\text{O} + \text{H}$



- Excellent agreement with previous work, but uncertainty in k substantially reduced

3.2 Methyl Ester Kinetics

1. Rate Constant Measurements

- OH + Methyl Formate
- OH + Methyl Acetate
- OH + Methyl Propanoate
- OH + Methyl Butanoate

Conducted as first-order experiments, with TBHP as source of OH



2. Methyl Formate Pyrolysis

- Multi-species Data

MF, CH₃OH
CO, CH₄, CH₂O

3.2 Methyl Ester Kinetics


- Summary: OH+Methyl Esters → Products



- Low scatter data with $\pm 25\%$ overall uncertainty
- Good agreement with modified SAR
(Structure Activity Relationship)

3.2 Methyl Ester Kinetics



- Comparison with Recent Quantum Calculations:
Methyl Formate + OH \rightarrow Products



- Shock tube/laser absorption rate 2-4x faster than calculation
- Confirms continuing value of high-accuracy experimental data

3.2 Methyl Ester Kinetics

- Methyl Formate Pyrolysis Kinetics:**
Multi-Species Laser Absorption Data Can Provide Near-Complete Oxygen Balance



@ $t = 300\mu\text{s}$



% Oxygen balance:

5.5% in CH_3OCHO

34.8% in CH_3OH

44.9% in CO

5.8% in CO_2

7.2% in CH_2O

Total: >98%

- Laser data successfully tracks all major contributors to O-atom balance
- Atom balance provides important new validation tool for chemical models

3.3 Butanol Isomer Kinetics


• Two Studies of Butanol Kinetics

1. n-Butanol Pyrolysis:

Challenge: Complex Pathways

Solution: Multi-Species Time-Histories


OH, CO, CH₄, H₂O, C₂H₄



2. tert-Butanol + OH → Products

➔ **Challenge:** Multiple Reaction Sites,
OH Reformation

Solution: Isotopic Labeling



3.3 Butanol Isomer Kinetics

- **OH⁺ tert-Butanol → Products**
Use of Isotopic Labeling
- **Challenges:** Multiple Reaction Sites, OH Reformation
- **Solution:** Isotopic labeling provides strategy to identify OH attack sites

3.3 Butanol Isomer Kinetics


- OH+ tert-Butanol → Products
Conventional Experiment using t-C₄H₉¹⁶OH

- Measurement of OH removal by t-C₄H₉¹⁶OH affected by OH reformation pathway

- Measurement gives

$$k_{1b} + k_{1a} (1 - BR_2)$$

$$= (k_{1b} + k_{1a}) (1 - BR_1 BR_2)$$



3.3 Butanol Isomer Kinetics

- OH+ tert-Butanol → Products
Conventional Experiment using t-C₄H₉¹⁶OH


- Measurement of OH removal by t-C₄H₉¹⁶OH affected by OH reformation pathway

- Measurement gives

$$k_{1b} + k_{1a}(1 - BR_2)$$


$$= (k_{1b} + k_{1a})(1 - BR_1 BR_2)$$

- Measurement of ¹⁶OH in t-C₄H₉¹⁸OH gives (k_{1a} + k_{1b}) with no ¹⁶OH reformation



3.3 Butanol Isomer Kinetics

- OH Absorption Data for tert-Butanol (O^{16} & O^{18})




- Slow removal of OH during ^{16}O butanol pyrolysis
- Faster removal of OH during ^{18}O butanol pyrolysis
- How do the rates compare?

3.3 Butanol Isomer Kinetics

- OH + tert-Butanol (O^{16} & O^{18})

Net OH Decay Rate due to Reaction with t-Butanol




$$k_{\text{total}} = (k_{1a} + k_{1b})$$

Derived from $t\text{-C}_4\text{H}_9^{18}\text{OH}$

3.3 Butanol Isomer Kinetics

- OH + tert-Butanol (O^{16} & O^{18})

Net OH Decay Rate due to Reaction with t-Butanol



$k_{\text{total}} = (k_{1a}+k_{1b})$
Derived from t-C₄H₉¹⁸OH

$(k_{1a}+k_{1b}) (1-BR_1BR_2)$
Derived from t-C₄H₉¹⁶OH

- First direct determination of overall rate ($k_{1a}+k_{1b}$)
- Ratio ($^{18}\text{O}/^{16}\text{O}$ expts.) gives $(1-BR_1BR_2) = 0.2$
- Since $BR_1 = 0.95$, then $BR_2 = 0.8 (\pm 0.05)$

4. New Species Diagnostics



Recent Work:

1. Aldehydes, e.g.: CH_2O , CH_3CHO
2. Alkenes, e.g.: iso- C_4H_8
3. Alkynes, e.g.: C_2H_2
4. Cavity Enhanced Absorption Spectroscopy (CEAS)

4. New Species Diagnostics - 1




- Aldehyde Diagnostics

Motivation: Aldehydes provide critical information about first stages of hydrocarbon oxidation, especially oxygenated fuels



Challenge: Overlapping absorption spectra


Strategy: Three colors (λ_1 and λ_2 in IR, and λ_3 in UV)




4. New Species Diagnostics - 1

- Validation Experiment: Simultaneous Measurement of Known CH_2O & CH_3CHO Concentrations
- Shock heat trioxane/acetaldehyde mixture
- Observe formaldehyde formation
- Recover correct $\text{CH}_2\text{O}/\text{CH}_3\text{CHO}$ concentrations

Measured Absorbance Time-Histories




Aldehyde Concentration Time-Histories



- Recovered $d\text{CH}_2\text{O}/dt$ & plateau mole fractions: Success!

4. New Species Diagnostics - 2

New Alkene Sensors: Far-IR Detection using QC Laser




- This spectral region holds high promise for multiple species: ethylene, propene, isobutene, ...

Example: isobutene, an important decomposition product of branched alkanes

4. New Species Diagnostics - 2


Iso-Octane Pyrolysis: First Detection of iso-Butene in Shock Tube



- Data provide important test for iso-octane decomposition pathways

4. New Species Diagnostics - 2

Iso-Octane Pyrolysis: Comparison with Modeling




- Data provide important test for iso-octane decomposition pathways
- Reasonable agreement with LLNL model!

Next alkyne detection: acetylene

4. New Species Diagnostics - 3

New Alkyne Sensor: Acetylene Detection at 3 μm




- **3 μm acetylene detection provides sensitive mid-IR detection**

Next: first application in propene

4. New Species Diagnostics - 3


First Acetylene Measurements: Propene Pyrolysis



- **Data provide high-sensitivity acetylene time-histories for modeling**

4. New Species Diagnostics - 3


First Acetylene Measurements: Propene Pyrolysis



- **Data provide high-sensitivity acetylene time-histories for modeling**

4. New Species Diagnostics - 3

First Acetylene Measurements: Propene Pyrolysis




- **Data provide high-sensitivity acetylene formation rates**
- **Reveals significant difference between JetSurF, ARAMCO and data**
- **Wide applicability of diagnostic: soot formation, alkene pyrolysis. ...**

Next: Cavity Enhanced Absorption Spectroscopy (CEAS)

4. New Species Diagnostics - 4

Ultra-sensitive Species Detection: Cavity-Enhanced Absorption Spectroscopy (CEAS) with TDL Source

➤ **Off-axis CEAS superior to traditional multipass**




- Cavity increases effective absorption path length, L_{eff}
- Variable gain via L_{eff} where $G=L_{\text{eff}}/D=R/(1-R)$ ($10 < G < 1000$ at low finesse)
- Off-axis insertion enables easy, robust alignment

Example: Detection limits for CO using MIR source

4. New Species Diagnostics - 4

CO formation by acetone (10ppm) pyrolysis



- **First highly-diluted CO formation experiment**
- **Detection limit $\sim 0.3\text{ppm}$ at 100kHz time resolution**
- **Data yield rate constant for acetone dissociation**

- **Demonstrates high CO sensitivity and good temporal resolution**
- **Further improvements: 10 – 100x in progress**

Summary:

Laser absorption & shock tubes are a frontier for combustion kinetics

Next: Three Lectures on Laser-Induced Fluorescence (LIF)

- ❖ Two level model
- ❖ More complex models
- ❖ Applications